

AIRL

BIO PHARMA

THREE STERILITY TESTING OPTIONS, ONE STANDARD OF EXCELLENCE

SCANRDI RAPID

Turnaround Time: 1 - 3 days

- Uses direct fluorescence and solid-phase laser scanning cytometry to detect microorganisms
- Validated alternative method per <1223> that is non-inferior to <71> testing
- Membrane Filtration
- Method Suitability Required (3.5x Sterility Volume)
- Can test filterable products only

CELSIS RAPID

Turnaround Time: 6 days

- Uses ATP-based bioluminescence to detect microorganisms
- Validated alternative method per <1223> that is non-inferior to <71> testing
- Closed Membrane Filtration or Direct Inoculation
- Method Suitability Required (4x Sterility Volume)
- Can test all sample types

USP <71>

Turnaround Time: 14 - 18 days

- Visual examination of growth media for turbidity
- Closed Membrane Filtration or Direct Inoculation
- Method Suitability Required (3x Sterility Volume)
- Can test all sample types